



silverglidesutures.com

Contents

SUTURES NON-ABSORBABLE

SHORTY™ 3 Fluoro Pink Polypropylene 4 Nylon 5 Fluoromid™

6

SUTURES ABSORBABLE

Polydioxanone Polyglycolic Acid Polyglycolic Acid Fast Absorbing 9 (Dental) **Polyglactin** 9 Glycolide & eCaprolactone 10

BOVINE SURGERY SUTURES

Fluoromid™ 11 Glycolide & eCaprolactone 11 (Uterus Closure) Polyglycolic Acid **Polydioxanone** 12 **Chromic Catgut Sutures** 12

NEEDLE CHART 13-14

SUTURE MATERIAL PROFILES 15





Dear Friends,

A warm welcome to our 2020 SilverGlide Veterinary Suture catalog.

September 2020 will mark our 6th anniversary.

During this period, our veterinary surgeons have been most forthcoming with regard to their daily frustrations, having to make do with sutures that have been designed primarily for human surgery. The insights they have provided have been instrumental in the re-design of traditional surgical sutures.

The result has been that during this 6 year period, SilverGlide has produced more than 60 sutures designed specifically for veterinary surgery; this works out to be a new suture release per month for each month that we have been in operation.

Most recent, and in response to requests for short sutures for securing wound drains and cannulas, or just to close a small wound without having to open a foil that contains a suture 70cm to 90cm long, we had great delight in launching SHORTY™ a LOW COST, HIGHLY VISIBLE new suture material (PVDF) that will knot securely with just 3 throws.

Further great ideas have been put to us and we can assure all our veterinary surgeons that more exciting developments are in the pipeline!

The Teasdel Family.

SHORTY™ Sutures

PV		POLYVINYLIDENE FLUORIDE, NON-ABSORBABLE							
ABSORPTION	N PROFILI	E	Exceptional Knot Security, knots flatten and lock against themselves, unwetted by blood, inert Non-Absorbable						
ORDER CODE	USP GAUGE	EP	NEEDLE TYPE	NEEDLE SIZE	NEEDLE	SUTURE LENGTH	QTY/ BOX	PRICE	
75514PVDF-2	3-0	2	3/8 Reverse Cutting $ abla$	25mm		45cm	36	\$81	
75514PVDF-1	2-0	3	3/8 Reverse Cutting $ abla$	25mm		45cm	36	\$81	

Double forward throw

(Surgeons Throw)



Single reverse throw

Creating a square knot, with which to adjust tension.



Single forward throw

(Locking the knot)



Final result

9 interrupted sutures, with long tails, from a 45cm length.



Attention All Veterinarians . . .



SHORTY

Sutures released!

Ever just needed a short length high vis suture for a small incision?

Look no further than SHORTY™.

Low Cost - No Waste.

Special introductory offer, for a limited time!

Innovation: Vets Talk, SilverGlide Listens


Fluoro Pink Polypropylene Sutures

RR			FI	LUORE	SCENT PINK MONOFILAMENT, NO	N-AB	SORB	ABLE
ABSORPTIO	N PROFILI	E	Unwetted by blood, inert, knots fl	atten and	lock against themselves for improved knot securit	y non-ae	SORBA	BLE
ORDER CODE	USP GAUGE	EP	NEEDLE TYPE	NEEDLE SIZE	NEEDLE	SUTURE LENGTH	OTY/ BOX	PRICE
75501RR-16	2-0	3	1/2 Circle Taper Point ⊙	26.5mm	\	70cm	36	\$190
75501RR-5	0	3.5	1/2 Olicio Tapol I olini	20.011111		70cm	36	\$190
75501RR-21	5-0	1	3/8 Circle Reverse Cutting ▽	16mm		45cm	12	\$69
75501RR-20	4-0	1.5	3/8 Circle Reverse Cutting ▽	18mm		45cm	12	\$69
75501RR-19	4-0	1.5	3/0 Circle neverse Culling V	10111111		45cm	36	\$138
75501RR-18	3-0	2				70cm	12	\$190
75501RR-9	3-0	2	3/8 Circle Reverse Cutting ▽	24mm		70cm	36	\$190
75501RR-1	3-0	2				45cm	36	\$138
75501RR-10	2-0	3				70cm	36	\$190
75501RR-2	2-0	3	3/8 Circle Reverse Cutting 🗸	30mm		45cm	36	\$138
75501RR-22	0	3.5				45cm	36	\$138
75501RR-12	0	3.5	1/2 Circle Reverse Cutting ▽	37mm		100cm	36	\$190
75501RR-13	0	3.5				100cm	36	\$190
75501RR-14	1	4				100cm	36	\$190
75501RR-15	2	5	3/8 Circle Reverse Cutting	50		100cm	36	\$190
75501RR-3	0	3.5	∇	50mm		45cm	36	\$138
75501RR-4	1	4				45cm	36	\$138
75501RR-23	2	5	-			45cm	36	\$138
75501RR-17	3-0	2	Straight Cutting 7	20	70cm	36	\$190	
75501RR-11	2-0	3	Straight Cutting $ abla$	60mm		70cm	36	\$190

[†] Indicated for Caslick Surgery and small incisions

Examples of high visibility sutures making for easier suture placement and removal.







Nylon Sutures

N					BLACK MONOFILAMENT, NO	N-AB	SORB	ABLE
ABSORPTIO	N PROFIL	E	High strength retention Non-Abs	orbable				
ORDER CODE	USP GAUGE	EP	NEEDLE TYPE	NEEDLE SIZE	NEEDLE	SUTURE LENGTH	OTY/ BOX	PRICE
75509N-30	6-0	0.7	3/8 Circle Reverse Cutting ▽	11mm	\smile	45cm	12	\$59
75509N-29	5-0	1	3/8 Circle Reverse Cutting ▽	12mm		45cm	12	\$59
75509N-28	5-0	1	0/0.01 1 1 1 1 1 1 1 1 1	40	\ /	45cm	12	\$59
75509N-27	5-0	1	- 3/8 Circle Reverse Cutting ▽	16mm		45cm	36	\$144
75509N-23	4-0	1.5	0/0.01 D	40	\ /	45cm	12	\$59
75509N-24	4-0	1.5	- 3/8 Circle Reverse Cutting ▽	18mm		45cm	36	\$144
75509N-25	4-0	1.5				45cm	12	\$59
75509N-26	4-0	1.5				70cm	36	\$144
75509N-18	3-0	2	3/8 Circle Reverse Cutting ▽	0.4		45cm	12	\$59
75509N-19	3-0	2		24mm		45cm	36	\$144
75509N-21	3-0	2				70cm	12	\$59
75509N-22	3-0	2				70cm	36	\$144
75509N-10	2-0	3			,	70cm	36	\$144
75509N-11	0	3.5	3/8 Circle Reverse Cutting ▽	30mm		100cm	36	\$144
75509N-15	1	4				100cm	36	\$144
75509N-12	1	4	3/8 Circle Reverse Cutting ▽	50mm		100cm	36	\$144
75509N-20	3-0	2	Otroinht Outling V	CO		70cm	36	\$144
75509N-17	2-0	3	- Straight Cutting ▽	60mm		70cm	36	\$144
75509N-16	2	5	3/8 Circle Reverse Cutting ▽	77mm		100cm	36	\$144

Braided Polyester Sutures

TieBack Equine

BP			BRAIDED PO	LYEST	ER, GREEN SILICONE COATED, NO	N-AB	SORB	ABLE	
ABSORPTION PRO	OFILE .		High strength braid Non-Absorbable						
ORDER CODE	USP GAUGE	EP	NEEDLE TYPE	NEEDLE SIZE	NEEDLE	SUTURE LENGTH	QTY/ BOX	PRICE	
75511BP-1	5	7	1/2 Circle Reverse Cutting ▽	37mm		75cm	12	\$59	
75511BP-2	5	7	1/2 Circle Taper Point ⊙	37mm		75cm	12	\$59	

Fluoromid™ Sutures

SMF			FLUORESCENT Y	ELLOW	MULTIFILAMENT POLYAMIDE, NO	ON-AB	SORB	ABLE		
ABSORPTION	N PROFIL	E	Easy handling, ultra high visibility Non-Absorbable							
ORDER CODE	USP GAUGE	EP	NEEDLE TYPE	NEEDLE SIZE	NEEDLE	SUTURE LENGTH	QTY/ BOX	PRICE		
75507SMF-8	3-0	2	3/8 Circle Reverse Cutting ▽	24mm		70cm	36	\$190		
75507SMF-9	2-0	3	3/8 Circle Reverse Cutting ▽	30mm		70cm	36	\$190		
75507SMF-14	0	3.5		JOHIII		100cm	36	\$190		
75507SMF-10	0	3.5	1/2 Circle Reverse Cutting $ riangle$	37mm		100cm	36	\$190		
75507SMF-11	1	4	3/8 Circle Reverse Cutting ▽	50mm		100cm	36	\$190		
75507SMF-15	2-0	3	Straight Cutting ▽	60mm		70cm	36	\$190		
75507SMF-16	3-0	2		60mm		70cm	36	\$190		



Contrasting yellow suture closure on a brindle Staffordshire Bull Terrier

Polydioxanone Sutures

PD0					PURPLE MONOFILAM	ENT, ABS	SORB	ABL
ABSORPTIO	N PROFILI	E	50-80% retained tensile strengt	h at 42 da y	ys. Essentially complete between 6-9 month	S.		
ORDER CODE	USP GAUGE	EP	NEEDLE TYPE	NEEDLE SIZE	NEEDLE	SUTURE LENGTH	OTY/ BOX	PRI
75513P0-17	5-0	1	1/2 Circle Taper Point ⊙	16mm	\bigcup	70cm	36	\$1
75513P0-15	4-0	1.5	1/2 Circle Taper Point ⊙	22mm		70cm	36	\$1
75513P0-12	3-0	2				70cm	36	\$1
75513P0-4	2-0	3	1/2 Circle Taper Point •	26.5mm		70cm	36	\$1
	0	3.5				100cm	36	\$1
75513P0-5	2-0	3				70cm	36	\$1
75513P0-9	1	4	1/2 Circle Taper Point •	36mm	100cm	36	\$1	
75513P0-8	0	3.5				100cm	36	\$1
75513P0-19	1	4				130cm	12	\$5
75513P0-11	2	5	1/2 Circle Taper Point	65mm		130cm	12	\$9
75513P0-21	3	6				100cm	12	\$9
75513PO-16	4-0	1.5	3/8 Circle Reverse Cutting ▽	18mm		70cm	36	\$1
75513P0-13	3-0	2	3/8 Circle Reverse Cutting	24mm		70cm	36	\$1
75513P0-14	2-0	3	3/8 Circle Reverse Cutting	26mm		70cm	36	\$1
75513P0-22	2-0	3	1/2 Circle Reverse Cutting ▽	37mm		70cm	36	\$1
75513PO-30	2-0	3	3/8 Circle Reverse Cutting ▽	40mm		45cm	36	\$1
75513P0-1	0	3.5				100cm	36	\$1
75513P0-2	1	4	1/2 Circle Reverse Cutting	40mm		100cm	36	\$1
75513P0-3	2	5				100cm	36	\$1
75513P0-23 ^t	0	3.5				45cm	36	\$1
75513P0-24 ^t	1	4	3/8 Circle Reverse Cutting 5			45cm	36	\$7
75513P0-25	2	5		F0~~		45cm	36	\$1
75513P0-28	0	3.5		50mm		100cm	36	\$
75513PO-18	1	4				100cm	36	\$
75513P0-20	2	5				100cm	36	\$

[†] Indicated for Caslick Surgery and small incisions

Polyglycolic Acid Sutures

PGA					BRAIDED VIOI	ET, AB	SORB	ABL
ABSORPTIO	N PROFILI	E	65% retained tensile strength at	14 days.	Essentially complete between 3-6 months .			
ORDER CODE	USP GAUGE	EP	NEEDLE TYPE	NEEDLE SIZE	NEEDLE	SUTURE LENGTH	OTY/ BOX	PRIC
75503PG-33	4-0	1.5	1/0.01 T	00		70cm	36	\$14
75503PG-32	3-0	2	1/2 Circle Taper Point ⊙	22mm		70cm	36	\$14
75503PG-25	3-0	2				70cm	36	\$14
75503PG-24	2-0	3			\	70cm	36	\$14
75503PG-10	0	3.5	1/2 Circle Taper Point ⊙	26.5mm	26.5mm	70cm	36	\$1
75503PG-30	1	4				100cm	36	\$1
75503PG-11	0	3.5			\	100cm	36	\$1
75503PG-12	1	4	1/2 Circle Taper Point ⊙	36mm		100cm	36	\$1
75503PG-17	0	3.5			\ /	100cm	36	\$1
75503PG-18	1	4	1/2 Circle Taper Point ⊙	40mm		100cm	36	\$1
75503PG-20	1	4				130cm	12	\$7
75503PG-21	2	5				130cm	12	\$7
75503PG-22	3	6	1/2 Circle Taper Point ⊙	65mm		130cm	12	\$7
75503PG-16	3	6				100cm	12	\$7
75503PG-31	3	6	3/8 Circle Taper Point ⊙	90mm		130cm	12	\$8
75503PG-34	6-0	0.7	2/0 Cirolo Dovorgo Cutting V	10mm	\ /	45cm	12	\$5
75503PG-28	5-0	1	3/8 Circle Reverse Cutting ▽	12mm		45cm	12	\$5
75503PG-27	4-0	1.5	3/8 Circle Reverse Cutting ▽	18mm		70cm	36	\$1
75503PG-26	3-0	2	3/8 Circle Reverse Cutting ▽	24mm	\ /	70cm	36	\$1
75503PG-1	3-0	2	3/0 OII OIE NEVELSE OULLING V	Z4111111		45cm	36	\$1
75503PG-23	2-0	3	1/2 Circle Reverse Cutting ▽	26mm		70cm	36	\$1
75503PG-3	0	3.5				45cm	36	\$1
75503PG-2	2-0	3	3/8 Circle Reverse Cutting ▽	30mm		45cm	36	\$1
75503PG-29	0	3.5				100cm	36	\$1
75503PG-4	1	4		F0		45cm	36	\$1
75503PG-36	2	5	3/8 Circle Reverse Cutting ▽	50mm		45cm	36	\$1

 $[\]dagger$ Indicated for Caslick Surgery and small incisions

Polyglycolic Acid Fast Absorbing Sutures

Dental

PGF	//	///		///	BRAIDED, UNDYE	D, AB	SORB	ABLE
ABSORPTION	N PROFIL	E	48% retained tensile strength at	7 days. E	ssentially complete at 14 days.			
ORDER CODE	USP GAUGE	EP	NEEDLE TYPE	NEEDLE SIZE	NEEDLE	SUTURE LENGTH	QTY/ BOX	PRICE
75503PG-38	5-0	1	1/2 Circle Reverse Cutting ▽	13mm	\bigcup	45cm	12	\$56
75503PG-7	4-0	1.5	3/8 Circle Reverse Cutting ▽	16mm		45cm	12	<i>\$56</i>
75503PG-37	3-0	2	3/8 Circle Reverse Cutting ▽	18mm	\smile	45cm	12	\$56



Polyglactin Sutures

PGLA				BRAIDED VIOLI	OLET, ABSORBABLE			
ABSORPTION	N PROFIL	E	74% retained tensile strength at	14 days.	Essentially complete at 63 days.			
ORDER CODE	USP GAUGE	EP	NEEDLE TYPE	NEEDLE SIZE	NEEDLE	SUTURE LENGTH	QTY/ BOX	PRICE
75503PGLA-1	3	6	1/2 Circle Taper Point ⊙	65mm		90cm	12	\$56



Glycolide & eCaprolactone Sutures

Dental

PGCL		PURPLE MONOFILAMENT, ABSORBABLE							
ABSORPTIO		E	60% retained tensile strength at	1	ssentially complete between 2-4 months.				
ORDER CODE	USP GAUGE	EP	NEEDLE TYPE	NEEDLE SIZE	NEEDLE	SUTURE LENGTH	QTY/ BOX	PRICE	
75505EC-17	5-0	1	1/2 Circle Taper Point ⊙	11mm	\bigcup	45cm	12	\$59	
75505EC-12	4-0	1.5	1/2 Circle Taper Point •	16mm		45cm	12	\$59	
75505EC-24	5-0	1	1/2 Circle Taper Point ⊙	18mm	\	45cm	12	\$59	
75505EC-22	4-0	1.5	1/2 Gilole Tapel Foliit	10111111		70cm	36	\$169	
75505EC-25	4-0	1.5			70cm	36	\$169		
75505EC-20	3-0	2	1/2 Circle Taper Point ⊙	22mm		70cm	36	\$169	
75505EC-19	2-0	3				70cm	36	\$169	
75505EC-26	3-0	2	1/2 Circle Taper Point ⊙ 26.5mm	00.5	\]	70cm	36	\$169	
75505EC-27	2-0	3		20.3111111		70cm	36	\$169	
75505EC-28	0	3.5	4 /0 O'l. T	00	\ /	100cm	36	\$169	
75505EC-29	1	4	1/2 Circle Taper Point ⊙	36mm		100cm	36	\$169	
75505EC-10	1	4	1/2 Circle Taper Point ⊙	65mm		130cm	12	\$85	
75505EC-23	5-0	1	3/8 Circle Reverse Cutting $ abla$	16mm		45cm	12	\$59	
75505EC-6	4-0	1.5	3/8 Circle Reverse Cutting $ abla$	18mm		70cm	36	\$169	
75505EC-21	3-0	2	3/8 Circle Reverse Cutting ▽	24mm		70cm	36	\$169	
75505EC-18	2-0	3	3/8 Circle Reverse Cutting ▽	26mm		70cm	36	\$169	
75505EC-7	2-0	3	1/0 0: 1 0 0	0.7	\ 1	70cm	36	\$169	
75505EC-30	0	3.5	1/2 Circle Reverse Cutting	37mm		100cm	36	\$169	

^{*} Bovine uterus closure only Codes in bold are ideal for Dental use.

Bovine Surgery

Fluoromid™ Sutures

SMF		FLUORESCENT YELLOW, NON-ABSORBABLE								
ABSORPTION	I PROFIL	E	Easy handling, ultra high visibility Non-Absorbable							
ORDER CODE	USP GAUGE	EP	NEEDLE TYPE	NEEDLE SIZE	NEEDLE	SUTURE LENGTH	QTY/ BOX	PRICE		
75507SMF-12	3	6	Straight Cutting $ abla$	65mm		130cm	12	\$90		
75507SMF-13	3	6	1/2 Circle Conventional Cutting / Straight Shank (Ski) △	75mm		130cm	12	\$90		

Glycolide & eCaprolactone Sutures

Uterus Closure

PGCL		PURPLE MONOFILAMENT, ABSORBABLE									
ABSORPTION	PROFILI	E	60% retained tensile strength at 7 days. Essentially complete between 2-4 months.								
ORDER CODE	USP GAUGE	EP	NEEDLE TYPE	NEEDLE SIZE	NEEDLE	SUTURE LENGTH	QTY/ BOX	PRICE			
75505EC-10	1	4	1/2 Circle Taper Point ⊙	65mm		130cm	12	\$85			

Polyglycolic Acid Sutures

PGA		BRAIDED VIOLET, ABSORBABLE							
ABSORPTIO	N PROFIL	OFILE 65% retained tensile strength at 14 days. Essentially complete between 3-6 months.							
ORDER CODE	USP GAUGE	EP	NEEDLE TYPE	NEEDLE SIZE	NEEDLE	SUTURE LENGTH	OTY/ BOX	PRICE	
75503PG-20	1	4				130cm	12	\$79	
75503PG-21	2	5	1/2 Circle Taper Point ⊙	65mm		130cm	12	\$79	
75503PG-22	3	6				130cm	12	\$79	
75503PG-31	3	6	3/8 Circle Taper Point ⊙	90mm		130cm	12	\$79	



Bovine Surgery

Polydioxanone Sutures

PD0		PURPLE MONOFILAMENT, ABSORBABLE							
ABSORPTION	N PROFIL	E	50-80% retained tensile strength at 42 days. Essentially complete between 6-9 months.						
ORDER CODE	USP GAUGE	EP	NEEDLE TYPE	NEEDLE SIZE	NEEDLE	SUTURE LENGTH	OTY/ BOX	PRICE	
75513P0-11	2	5	1/2 Circle Taper Point ⊙	65mm		130cm	12	\$91	

Chromic Catgut Sutures

CG		NATURAL ABSORBABLE							
ABSORPTION	ION PROFILE 100% retained tensile strength at 7 days. Essentially complete in 3 months.								
ORDER CODE	USP GAUGE	EP	NEEDLE TYPE	NEEDLE SIZE	NEEDLE	SUTURE LENGTH	OTY/ BOX	PRICE	
75515CG-1	1	4	1/0 Cirolo Topor Doint	GEmm		130cm	12	\$56	
75515CG-2	2	5	1/2 Circle Taper Point ⊙	65mm		130cm	12	\$56	



Circle Taper Point & Straight Cutting

1/2 Circle Taper Point •

Needle	Size	SilverGlide™	Ethicon™ (similar)	Syneture™	B.Braun ™ (similar)
$\overline{}$	● 11mm	T11	-		
\bigcup_{\cdot}	● 13mm	T22	TF	CV-22	HR13
	● 16mm	T16	-		
	◆ 18mm	T28	RB-1	CV-23	HR18
	◆ 22mm	T4	SH-1	CV-25	HR22
	⊙ 26.5mm	T5	SH / CT-1	V20 / GS22	HR26
	⊙ 36mm	Т8	MH / V-26	V-26	HR36
	⊙ 37mm	T37	CT-1	-	-
	● 40mm	T40	-	-	-
	⊙ 65mm	Т39	TP-1	-	-

3/8 Circle Taper Point \odot



Straight Cutting ∇

Needle	Size	SilverGlide™	Ethicon ™ (similar)	Syneture™	B.Braun ™ (similar)
Straight Cutting	▽ 60mm	CE60	KS	SC-2	GS60
Straight Cutting	▽ 65mm	CE65	KS	SC-2	GS60
1/2 Circle Conventional Cutting / Straight Shank (Ski) △	<u></u> 75mm	SKI75	-	-	-

Circle Reverse Cutting

1/2 Circle Reverse Cutting ∇

Needle	Size	SilverGlide™	Ethicon™ (similar)	Syneture™	B.Braun ™ (similar)
	√ 13mm	C13	-	-	-
	▽ 26mm	C21	CP-2	GS-10	HS26
	▽ 37mm	C37	CP-1	GS-11	HS37
	▽ 40mm	C14	СР	GS-12	HS40

3/8 Circle Reverse Cutting ∇

Needle 	Size	SilverGlide™	Ethicon™ (similar)	Syneture™	B.Braun ™ (similar)
	▽11mm	C16	C-2	C-1	DS12
	▽12mm	C17	C-2	C-1	DS12
	▽16mm	C22	FS-3	C-12	DS16
	▽18mm	C6	FS-2	C-13	DS18
	▽ 24mm	C7	FS-1	C-14	DS24
	▽ 25mm	C25	-	-	-
	▽ 26mm	C26	FS	-	-
	▽30mm	C10	FSL	C-16	DS30
	▽ 40mm	C14	FSLX	C-17	
	▽ 50mm	C50	-	-	-
	▽ 77mm	C77	LR	-	-

Suture Material Absorption Profile

NON-ABSORBABL	E	MATERIAL DESCRIPTION	QU	ALITIES	SIMILAR TO	
NYLON		Black monofilament	High strength retention	High strength retention		
FLUOROMID™		Fluorescent Yellow polyamide multifilament	Easy handling, ultra high visibi Fluorescent yellow			
WHITE FLUOROMID™		White polyamide multifilament	Easy handling			
SHORTY™ POLYVINYLIDENE FLUORIDE YELLOW		Bright Yellow Monofilament	High visibilty Exceptional Knot Security, knothemselves, Unwetted by blood			
BLUE POLYPROPYLENE		Dark blue monofilament	Unwetted by blood, inert, knots for improved knot security	Prolene [™] * Surgipro [™] §/ Surgilene [™] § Premilene [™] †		
FLUORESCENT PINK POLYPROPYLENE		Fluorescent Pink monofilament	Ultra high visibility Fluorescent pink Unwetted by blood, inert, knots for improved knot security			
POLYESTER		Braided polyester green silicone coated	High strength braid	Ticron § Ethibond *		
ABSORBABLE		MATERIAL DESCRIPTION	BSR PROFILE	MASS ABSORPTION	SIMILAR TO	
GLYCOLIDE & ECAP		Monofilament purple	In vivo strength retention 60% at 7 days 25% at 14 days	Essentially complete between 2-4 months	Monocryl™ * Caprosyn™ §/ Biosyn™ § Monosyn™ †	
POLYGLYCOLIC ACID		Braided violet	In vivo strength retention 65% at 14 days 20% at 21 days	Essentially complete between 3-6 months	Vicryl™* Safil™† Polsorb§	
POLYGLYCOLIC ACID FAST ABSORBING		Braided undyed	, ,		Vicryl Rapide * Safil Quick †	
POLYGLACTIN		Braided violet	In vivo strength retention 74% at 14 days 9 weeks 49% at 21 days 24% at 28 days		Vicryl™	
POLYDIOXANONE		Monofilament purple	In vivo strength retention 80-90% at 14 days 60-80% at 28 days 50-80% at 42 days	Essentially complete between 6-9 months	PDS11 [™] * MonoPlus [™] † Maxon §	

^{*} Trademark of Ethicon § Trademark of Covidien † Trademark of B.Braun



For sample requests or to place an order call **800 447 0687** or visit our website **silverglidesutures.com**

Distributed by:



Phone: 800 447 0687 **Fax:** 1 715 953 2731

Email: sales@maianimalhealth.com

MAI Animal Health

Address: 605 Pro-Ject Drive, Elmwood, WI 54740