



**AUSTRALIAN  
VETERINARY  
CATALOGUE  
& PRICE LIST  
JANUARY 2020**

***[silverglide.com.au](http://silverglide.com.au)***

# Contents

## SUTURES NON-ABSORBABLE

<b>SHORTY™</b>	3
Fluoro Pink Polypropylene	4
Polypropylene	5
Braided Polyester	5
Nylon	6
Fluoromid™	7
White Fluoromid™ Cassettes	7

## SUTURES ABSORBABLE

Polydioxanone	8
Polyglycolic Acid	9
Polyglycolic Acid Fast Absorbing (Dental)	10
Polyglactin	10
Glycolide & eCaprolactone	11

## SUTURES NON-NEEDLE ABSORBABLE

Polyglycolic Acid	12
Glycolide & eCaprolactone	12
Polydioxanone	12
Reusable Stapler	12

## BOVINE SURGERY SUTURES

Fluoromid™	13
Glycolide & eCaprolactone (Uterus Closure)	13
Polyglycolic Acid	13
Polydioxanone	13

## SKIN CLOSURE

SurgiSeal Vet™ (2-Octyl Skin Adhesive)	14
--	----

## CORNEAL EYE DROPS

Clerapliq®	15
------------	----

## PAW FRICTION™

16

## VITAL SIGNS MONITOR

VetCorder	17
-----------	----

## PATIENT WARMING

Hotdog®	18
---------	----

## NEEDLE CHART

19-22

## SUTURE MATERIAL PROFILES

23



*Dear Friends,*

*A warm welcome to our 2020 SilverGlide Veterinary Suture catalogue.*

*September 2019 marked our 5th anniversary.*

*During this period, our veterinary surgeons have been most forthcoming with regard to their daily frustrations, having to make do with sutures that have been designed primarily for human surgery. The insights they have provided have been instrumental in the re-design of traditional surgical sutures.*

*The result has been that during this 5 year period, SilverGlide has produced more than 60 sutures designed specifically for veterinary surgery; this works out to be a new suture release per month for each month that we have been in operation.*

*Most recent, and in response to requests for short sutures for securing wound drains and cannulas, or just to close a small wound without having to open a foil that contains a suture 70cm to 90cm long, we had great delight in launching SHORTY™ a LOW COST, HIGHLY VISIBLE new suture material (PVDF) that will knot securely with just 3 throws.*

*Further great ideas have been put to us and we can assure all our veterinary surgeons that more exciting developments are in the pipeline!*

***The Teasdel Family.***



### Suture Code Rationale

Our suture codes are designed to describe exactly what they are...  
Let's take SHORTY™ as an example:

Suture order code :	PV	=	Polyvinylidene Fluoride suture material
<b>PV30425D</b>	<b>30</b>	=	3-0 usp gauge
	<b>4</b>	=	45cm suture length
	<b>25</b>	=	25mm needle
	<b>D</b>	=	36 sutures per box


# Sutures *Non-Absorbable*

## SHORTY™ Sutures

PV		POLYVINYLIDENE FLUORIDE, NON-ABSORBABLE							
ABSORPTION PROFILE			Exceptional Knot Security, knots flatten and lock against themselves, unwetted by blood, inert <b>Non-Absorbable</b>						
ORDER CODE	USP GAUGE	EP	NEEDLE TYPE	NEEDLE SIZE	NEEDLE	SUTURE LENGTH	QTY/ BOX	PRICE	
PV30425D	3-0	2	3/8 Reverse Cutting ▽	25mm		45cm	36	\$90	
PV20425D	2-0	3	3/8 Reverse Cutting ▽	25mm		45cm	36	\$90	


### Double forward throw

(Surgeons Throw)




### Single reverse throw

Creating a square knot, with which to adjust tension.




### Single forward throw

(Locking the knot)



### Final result

9 interrupted sutures, with long tails, from a 45cm length.



# Attention All Veterinarians . . .



SHORTY™ DEALS

## SHORTY™ Sutures released!

Ever just needed a short length high vis suture for a small incision?

**Look no further than SHORTY™.**

■ Low Cost - No Waste.









Special introductory offer, for a limited time!

Innovation: Vets Talk, SilverGlide Listens



# Sutures *Non-Absorbable*

## Fluoro Pink Polypropylene Sutures

RR		FLUORESCENT PINK MONOFILAMENT, NON-ABSORBABLE						
ABSORPTION PROFILE			Unwetted by blood, inert, knots flatten and lock against themselves for improved knot security <b>NON-ABSORBABLE</b>					
ORDER CODE	USP GAUGE	EP	NEEDLE TYPE	NEEDLE SIZE	NEEDLE	SUTURE LENGTH	QTY/ BOX	PRICE
RR2075D	2-0	3	1/2 Circle Taper Point ☉	26.5mm		70cm	36	\$154
RR0075D	0	3.5				70cm	36	\$154
RR50422N	5-0	1	3/8 Circle Reverse Cutting ▽	16mm		45cm	12	\$55
RR4046N	4-0	1.5	3/8 Circle Reverse Cutting ▽	18mm		45cm	12	\$55
RR4046D	4-0	1.5				45cm	36	\$154
RR3077N	3-0	2	3/8 Circle Reverse Cutting ▽	24mm		70cm	12	\$55
RR3077D	3-0	2				70cm	36	\$154
RRCLK3047D†	3-0	2				45cm	36	\$110
RR20710D	2-0	3	3/8 Circle Reverse Cutting ▽	30mm		70cm	36	\$154
RRCLK20410D†	2-0	3				45cm	36	\$110
RRCLK00310D	0	3.5				45cm	36	\$110
RR00113D	0	3.5	1/2 Circle Reverse Cutting ▽	37mm		100cm	36	\$154
RR00130D	0	3.5	3/8 Circle Reverse Cutting ▽	50mm		100cm	36	\$154
RR01130D	1	4				100cm	36	\$154
RR02130D	2	5				100cm	36	\$154
RRCLK00430D†	0	3.5				45cm	36	\$110
RRCLK01430D†	1	4				45cm	36	\$110
RRCLK02430D†	2	5				45cm	36	\$110
RR30760D	3-0	2	Straight Cutting ▽	60mm		70cm	36	\$154
RR20760D	2-0	3				70cm	36	\$154



† Indicated for Caslick Surgery and small incisions

***Examples of high visibility sutures making for easier suture placement and removal.***





# Sutures *Non-Absorbable*


## Polypropylene Sutures

PP		BLUE MONOFILAMENT, NON-ABSORBABLE							
ABSORPTION PROFILE		Unwetted by blood, inert, knots flatten and lock against themselves for improved knot security <b>Non-Absorbable</b>							
ORDER CODE	USP GAUGE	EP	NEEDLE TYPE	NEEDLE SIZE	NEEDLE	SUTURE LENGTH	QTY/ BOX	PRICE	
PP3077D	3-0	2	3/8 Circle Reverse Cutting ▽	24mm		70cm	36	\$154	
PP20710D	2-0	3	3/8 Circle Reverse Cutting ▽	30mm		70cm	36	\$154	

## Braided Polyester Sutures










*TieBack Equine*

BP		BRAIDED POLYESTER, GREEN SILICONE COATED, NON-ABSORBABLE							
ABSORPTION PROFILE		High strength braid <b>Non-Absorbable</b>							
ORDER CODE	USP GAUGE	EP	NEEDLE TYPE	NEEDLE SIZE	NEEDLE	SUTURE LENGTH	QTY/ BOX	PRICE	
BP05737N	5	7	1/2 Circle Reverse Cutting ▽	37mm		75cm	12	\$88	
BP05738N	5	7	1/2 Circle Taper Point ⊙	37mm		75cm	12	\$88	








# Sutures *Non-Absorbable*

## Nylon Sutures

N			BLACK MONOFILAMENT, NON-ABSORBABLE						
ABSORPTION PROFILE			High strength retention <b>Non-Absorbable</b>						
ORDER CODE	USP GAUGE	EP	NEEDLE TYPE	NEEDLE SIZE	NEEDLE	SUTURE LENGTH	QTY/ BOX	PRICE	
N60416N	6-0	0.7	3/8 Circle Reverse Cutting ▽	11mm		45cm	12	\$44	
N50417N	5-0	1	3/8 Circle Reverse Cutting ▽	12mm		45cm	12	\$44	
N50422N	5-0	1	3/8 Circle Reverse Cutting ▽	16mm		45cm	12	\$44	
N50422D	5-0	1				45cm	36	\$110	
N4046N	4-0	1.5	3/8 Circle Reverse Cutting ▽	18mm		45cm	12	\$44	
N4046D	4-0	1.5				45cm	36	\$110	
N4047N	4-0	1.5	3/8 Circle Reverse Cutting ▽	24mm		45cm	12	\$44	
N4077D	4-0	1.5				70cm	36	\$110	
N3047N	3-0	2				45cm	12	\$44	
N3047D	3-0	2				45cm	36	\$110	
N3077D	3-0	2				70cm	36	\$110	
N3077N	3-0	2				70cm	12	\$44	
N20710D	2-0	3				70cm	36	\$110	
N00110D	0	3.5	3/8 Circle Reverse Cutting ▽	30mm		100cm	36	\$110	
N01140D	1	4				100cm	36	\$110	
N01130D	1	4	3/8 Circle Reverse Cutting ▽	50mm		100cm	36	\$110	
N30760D	3-0	2	Straight Cutting ▽	60mm		70cm	36	\$110	
N20760D	2-0	3				70cm	36	\$110	
N02125D	2	5	3/8 Circle Reverse Cutting ▽	77mm		100cm	36	\$154	

# Sutures *Non-Absorbable*

## Fluoromid™ Sutures

SMF		FLUORESCENT YELLOW MULTIFILAMENT POLYAMIDE, NON-ABSORBABLE							
ABSORPTION PROFILE			Easy handling, ultra high visibility <b>Non-Absorbable</b>						
ORDER CODE	USP GAUGE	EP	NEEDLE TYPE	NEEDLE SIZE	NEEDLE	SUTURE LENGTH	QTY/ BOX	PRICE	
SMF3077D	3-0	2	3/8 Circle Reverse Cutting ▽	24mm		70cm	36	\$110	
SMF20710D	2-0	3	3/8 Circle Reverse Cutting ▽	30mm		70cm	36	\$110	
SMF00110D	0	3.5				100cm	36	\$110	
SMF00113D	0	3.5	1/2 Circle Reverse Cutting ▽	37mm		100cm	36	\$110	
SMF01130D	1	4	3/8 Circle Reverse Cutting ▽	50mm		100cm	36	\$110	
SMF20760D	2-0	3	Straight Cutting ▽	60mm		70cm	36	\$110	
SMF30760D	3-0	2				70cm	36	\$110	













## White Fluoromid™ Cassette

White Fluoromid™				
ABSORPTION PROFILE		Easy handling <b>Non-Absorbable</b>		
CODE	USP GAUGE	LENGTH	QTY	PRICE
SMCAS40100	4-0	100	1	\$66
SMCAS30100	3-0	100	1	\$66
SMCAS20100	2-0	100	1	\$66
SMCAS00100	0	100	1	\$66
SMCAS01050	1	50	1	\$66
SMCAS02050	2	50	1	\$66



# Sutures *Absorbable*

## Polydioxanone Sutures

























PDO			PURPLE MONOFILAMENT, ABSORBABLE						
ABSORPTION PROFILE			50-80% retained tensile strength at <b>42 days</b> . Essentially complete between <b>6-9 months</b> .						
ORDER CODE	USP GAUGE	EP	NEEDLE TYPE	NEEDLE SIZE	NEEDLE	SUTURE LENGTH	QTY/ BOX	PRICE	
P050724D	5-0	1	1/2 Circle Taper Point ◉	16mm		70cm	36	\$194	
P04074D	4-0	1.5	1/2 Circle Taper Point ◉	22mm		70cm	36	\$194	
P03075D	3-0	2	1/2 Circle Taper Point ◉	26.5mm		70cm	36	\$194	
P02075D	2-0	3				70cm	36	\$194	
P00015D	0	3.5				100cm	36	\$194	
P02078D	2-0	3	1/2 Circle Taper Point ◉	36mm		70cm	36	\$194	
P00118D	1	4				100cm	36	\$194	
P00018D	0	3.5				100cm	36	\$194	
P0011339N	1	4	1/2 Circle Taper Point ◉	65mm		130cm	12	\$80	
P0021339N	2	5				130cm	12	\$80	
P003139N	3	6				100cm	12	\$80	
P04076D	4-0	1.5	3/8 Circle Reverse Cutting ▽	18mm		70cm	36	\$194	
P03077D	3-0	2	3/8 Circle Reverse Cutting ▽	24mm		70cm	36	\$194	
P020726D	2-0	3	3/8 Circle Reverse Cutting ▽	26mm		70cm	36	\$194	
P020713D	2-0	3	1/2 Circle Reverse Cutting ▽	37mm		70cm	36	\$194	
POCLK20410D†	2-0	3	3/8 Circle Reverse Cutting ▽	40mm		45cm	36	\$134	
P000114D	0	3.5	1/2 Circle Reverse Cutting ▽	40mm		100cm	36	\$194	
P001114D	1	4				100cm	36	\$194	
P002114D	2	5				100cm	36	\$194	
P001130D	1	4	3/8 Circle Reverse Cutting ▽	50mm		100cm	36	\$194	
P002130D	2	5				100cm	36	\$194	
POCLK00430D†	0	3.5				45cm	36	\$134	
POCLK01430D†	1	4				45cm	36	\$134	
POCLK02430D†	2	5				45cm	36	\$134	
P0001130D	0	3.5				100cm	36	\$194	

† Indicated for Caslick Surgery and small incisions



# Sutures *Absorbable*

## Polyglycolic Acid Sutures




PGA			BRAIDED VIOLET, ABSORBABLE					
ABSORPTION PROFILE			65% retained tensile strength at 14 days. Essentially complete between 3-6 months.					
ORDER CODE	USP GAUGE	EP	NEEDLE TYPE	NEEDLE SIZE	NEEDLE	SUTURE LENGTH	QTY/ BOX	PRICE
PG4074D	4-0	1.5	1/2 Circle Taper Point 	22mm		70cm	36	\$174
PG3074D	3-0	2				70cm	36	\$174
PG3075D	3-0	2	1/2 Circle Taper Point 	26.5mm		70cm	36	\$174
PG2075D	2-0	3				70cm	36	\$174
PG0075D	0	3.5				70cm	36	\$174
PG0115D	1	4				100cm	36	\$174
PG0018D	0	3.5	1/2 Circle Taper Point 	36mm		100cm	36	\$174
PG0118D	1	4				100cm	36	\$174
PG00140D	0	3.5	1/2 Circle Taper Point 	40mm		100cm	36	\$174
PG01140D	1	4				100cm	36	\$174
PG011339N	1	4	1/2 Circle Taper Point 	65mm		130cm	12	\$80
PG021339N	2	5				130cm	12	\$80
PG031339N	3	6				130cm	12	\$80
PG03139N	3	6				100cm	12	\$80
PG031390N	3	6	3/8 Circle Taper Point 	90mm		130cm	12	\$90
PG60417N	6-0	0.7	3/8 Circle Reverse Cutting 	12mm		45cm	12	\$96
PG50417N	5-0	1				45cm	12	\$96
PG4076D	4-0	1.5	3/8 Circle Reverse Cutting 	18mm		70cm	36	\$174
PG3077D	3-0	2	3/8 Circle Reverse Cutting 	24mm		70cm	36	\$174
PGCLK3047D†	3-0	2				45cm	36	\$137
PG20721D	2-0	3	1/2 Circle Reverse Cutting 	26mm		70cm	36	\$174
PG00110D	0	3.5	3/8 Circle Reverse Cutting 	30mm		100cm	36	\$174
PGCLK20410D†	2-0	3				45cm	36	\$134
PGCLK00410D†	0	3.5				45cm	36	\$134
PGCLK01430D†	1	4	3/8 Circle Reverse Cutting 	50mm		45cm	36	\$134
PGCLK02430D†	2	5				45cm	36	\$134

† Indicated for Caslick Surgery and small incisions

# Sutures *Absorbable*


## Polyglycolic Acid Fast Absorbing Sutures

*Dental*

PGF		BRAIDED, UNDYED, ABSORBABLE							
ABSORPTION PROFILE			48% retained tensile strength at <b>7 days</b> . Essentially complete at <b>14 days</b> .						
ORDER CODE	USP GAUGE	EP	NEEDLE TYPE	NEEDLE SIZE	NEEDLE	SUTURE LENGTH	QTY/ BOX	PRICE	
PGF50413N	5-0	1	1/2 Circle Reverse Cutting ▽	13mm		45cm	12	\$77	
PGF40422N	4-0	1.5	3/8 Circle Reverse Cutting ▽	16mm		45cm	12	\$77	
PGF3046N	3-0	2	3/8 Circle Reverse Cutting ▽	18mm		45cm	12	\$77	















## Polyglactin Sutures

PGLA		BRAIDED VIOLET, ABSORBABLE							
ABSORPTION PROFILE			74% retained tensile strength at <b>14 days</b> . Essentially complete at <b>63 days</b> .						
ORDER CODE	USP GAUGE	EP	NEEDLE TYPE	NEEDLE SIZE	NEEDLE	SUTURE LENGTH	QTY/ BOX	PRICE	
PGLA03939N	3	6	1/2 Circle Taper Point ⊙	65mm		90cm	12	\$99	



# Sutures *Absorbable*

## Glycolide & eCaprolactone Sutures

GCL			PURPLE MONOFILAMENT, ABSORBABLE						
ABSORPTION PROFILE			60% retained tensile strength at 7 days. Essentially complete between 2-4 months.						
ORDER CODE	USP GAUGE	EP	NEEDLE TYPE	NEEDLE SIZE	NEEDLE	SUTURE LENGTH	QTY/ BOX	PRICE	
EC50450N	5-0	1	1/2 Circle Taper Point ◉	11mm		45cm	12	\$86	
EC40416N	4-0	1.5	1/2 Circle Taper Point ◉	16mm		45cm	12	\$86	
EC50428N	5-0	1	1/2 Circle Taper Point ◉	18mm		45cm	12	\$86	
EC40728D	4-0	1.5				70cm	36	\$184	
EC4074D	4-0	1.5	1/2 Circle Taper Point ◉	22mm		70cm	36	\$184	
EC3074D	3-0	2				70cm	36	\$184	
EC2074D	2-0	3				70cm	36	\$184	
EC3075D	3-0	2	1/2 Circle Taper Point ◉	26.5mm		70cm	36	\$184	
EC2075D	2-0	3				70cm	36	\$184	
EC0018D	0	3.5	1/2 Circle Taper Point ◉	36mm		100cm	36	\$184	
EC0118D	1	4				100cm	36	\$184	
EC011339N*	1	4	1/2 Circle Taper Point ◉	65mm		130cm	12	\$80	
EC50416N	5-0	1	3/8 Circle Reverse Cutting ▽	16mm		45cm	12	\$86	
EC4076D	4-0	1.5	3/8 Circle Reverse Cutting ▽	18mm		70cm	36	\$184	
EC3077D	3-0	2	3/8 Circle Reverse Cutting ▽	24mm		70cm	36	\$184	
EC20726D	2-0	3	3/8 Circle Reverse Cutting ▽	26mm		70cm	36	\$184	
EC20713D	2-0	3	1/2 Circle Reverse Cutting ▽	37mm		70cm	36	\$184	
EC00113D	0	3.5				100cm	36	\$184	

\* Bovine uterus closure only

# Sutures Non-Needled *Absorbable*


Polyglycolic Acid				
ABSORPTION PROFILE	65% retained tensile strength at 14 days. Essentially complete between 3-6 months			
CODE	USP GAUGE	LENGTH	QTY/ BOX	PRICE
PGNN3045D	3-0	45	36	\$79
PGNN2045D	2-0	45	36	\$79
PGNN0045D	0	45	36	\$79
PGNN0145D	1	45	36	\$79
PGNN0245D	2	45	36	\$79
PGNN30100D	3-0	100	36	\$119
PGNN20100D	2-0	100	36	\$119
PGNN00100D	0	100	36	\$119
PGNN01100D	1	100	36	\$132
PGNN02100D	2	100	36	\$132

Polydioxanone				
ABSORPTION PROFILE	50-80% retained tensile strength at 42 days. Essentially complete between 6-9 months			
CODE	USP GAUGE	LENGTH	QTY/ BOX	PRICE
PONN3045D	3-0	45	36	\$79
PONN2045D	2-0	45	36	\$79
PONN0045D	0	45	36	\$79
PONN0145D	1	45	36	\$79
PONN0245D	2	45	36	\$79
PONN30100D	3-0	100	36	\$119
PONN20100D	2-0	100	36	\$119
PONN00100D	0	100	36	\$119
PONN01100D	1	100	36	\$132
PONN02100D	2	100	36	\$132

Glycolide & eCap				
ABSORPTION PROFILE	60% retained tensile strength at 7 days Essentially complete between 2-4 months			
CODE	USP GAUGE	LENGTH	QTY/ BOX	PRICE
ECNN3045D	3-0	45	36	\$79
ECNN2045D	2-0	45	36	\$79
ECNN0045D	0	45	36	\$79
ECNN0145D	1	45	36	\$79
ECNN0245D	2	45	36	\$79
ECNN30100D	3-0	100	36	\$119
ECNN20100D	2-0	100	36	\$119
ECNN00100D	0	100	36	\$119
ECNN01100D	1	100	36	\$132
ECNN02100D	2	100	36	\$132



## Reusable Staplers

CODE	DESCRIPTION	QUANTITY	PRICE
SRH-1	Pre-Loaded SRH Handle	1	\$36
SRC35W	Cartridges	1	\$19.50
SRE	Staple Extractor	1	\$14.50




# Bovine Surgery

## Fluoromid™ Sutures



SMF		FLUORESCENT YELLOW OR WHITE, NON-ABSORBABLE							
ABSORPTION PROFILE			Easy handling, ultra high visibility <b>Non-Absorbable</b>						
ORDER CODE	USP GAUGE	EP	NEEDLE TYPE	NEEDLE SIZE	NEEDLE	SUTURE LENGTH	QTY/ BOX	PRICE	
SMF031360N	3	6	Straight Cutting ▽	65mm		130cm	12	\$66	
SMF031375N	3	6	1/2 Circle Conventional Cutting / Straight Shank (Ski) △	75mm		130cm	12	\$66	

## Glycolide & eCaprolactone Sutures


Uterus Closure

GCL		PURPLE MONOFILAMENT, ABSORBABLE							
ABSORPTION PROFILE			<b>60%</b> retained tensile strength at <b>7 days</b> . Essentially complete between <b>2-4 months</b> .						
ORDER CODE	USP GAUGE	EP	NEEDLE TYPE	NEEDLE SIZE	NEEDLE	SUTURE LENGTH	QTY/ BOX	PRICE	
EC011339N	1	4	1/2 Circle Taper Point ⊙	65mm		130cm	12	\$80	

## Polyglycolic Acid Sutures

PGA		BRAIDED VIOLET, ABSORBABLE							
ABSORPTION PROFILE			<b>65%</b> retained tensile strength at <b>14 days</b> . Essentially complete between <b>3-6 months</b> .						
ORDER CODE	USP GAUGE	EP	NEEDLE TYPE	NEEDLE SIZE	NEEDLE	SUTURE LENGTH	QTY/ BOX	PRICE	
PG011339N	1	4	1/2 Circle Taper Point ⊙	65mm		130cm	12	\$80	
PG021339N	2	5	1/2 Circle Taper Point ⊙	65mm		130cm	12	\$80	
PG031339N	3	6	1/2 Circle Taper Point ⊙	65mm		130cm	12	\$80	
PG031390N	3	6	3/8 Circle Taper Point ⊙	90mm		130cm	12	\$90	

## Polydioxanone Sutures

PDO		PURPLE MONOFILAMENT, ABSORBABLE							
ABSORPTION PROFILE			<b>50-80%</b> retained tensile strength at <b>42 days</b> . Essentially complete between <b>6-9 months</b> .						
ORDER CODE	USP GAUGE	EP	NEEDLE TYPE	NEEDLE SIZE	NEEDLE	SUTURE LENGTH	QTY/ BOX	PRICE	
P0021339N	2	5	1/2 Circle Taper Point ⊙	65mm		130cm	12	\$80	



# Skin Closure


## SurgiSeal Vet™ 2-Octyl Skin Adhesive

2-Octyl Cyanoacrylate. Flexible, Gentle, Waterproof, Pathogen Barrier			
CODE	DESCRIPTION	QUANTITY	PRICE
SV0350B	3.5ml. 1 bottle, 10 foam applicator swabs, 2 pipettes (applicator pack)	1	\$39
SV012A	20 foam applicator swabs, 4 pipettes (applicator pack)	1	\$11
SS050S	0.50ml Single Use applicators <b>NEW LARGER APPLICATOR</b>	6	\$72

- Layered along and lateral to the wound creating wound edge support
- SurgiSeal Vet™ is as flexible as skin; No cracking
- Immediately waterproof
- Water vapour & oxygen transmissible
- Traps and kills skin bacteria, preventing pathogen migration to the wound
- Protects suture or staple lines from post-operative contamination
- ECONOMICAL – Just \$1-2 per application

### COMPARISONS BETWEEN 2-OCTYL CYANACROLATE & BUTYL CYANACROLATE


#### Greater Breaking Strength<sup>1,2</sup>



[1]Perry LC. An evaluation of acute incisional strength with traumaseal surgical tissue adhesive wound closure. Findings by: Dimensional Analysis Systems, Inc., Leonia, NJ. 1995.  
 [2] Singer AJ, et al. Comparison of wound-bursting strengths and surface characteristics of FDA-approved tissue adhesive for skin closure. J Adhes Sci Technol. 2004;18:19-27.

### APPLICATION DIFFERENCES BETWEEN 2-OCTYL CYANACROLATE & BUTYL CYANACROLATE

**Octyl layering** - Is layered on and has a wide area of application. This **wider bridging** action also **contributes** to its **increased wound-breaking strength** and **antimicrobial properties** by immobilizing skin flora.



**Butyl beading** - The butyl cyanoacrylates form a thin, narrow bridge over the wound compared with octyl cyanoacrylates




# Corneal Eye Drops

## Clerapliq®

Corneal Ulcer Drops. Enables Regrowth of Epithelium & Significantly Reduces Pain			
CODE	DESCRIPTION	QUANTITY	PRICE
CPLQ1	Clerapliq box of 5 pipette doses	1	\$69


**Clerapliq® breakthrough technology for the treatment of pet's persistent or chronic corneal wounds.**



5 year old Jack Russel.  
Diagnosis epithelial ulcer using fluoroscene.



ONE DROP OF CLERAPLIQ  
EVERY SECOND DAY.



Complete healing.

CLERAPLIQ® is an ophthalmologic matrix therapy agent. It triggers regeneration of damaged tissue by restoring the corneal matrix architecture, thus encouraging corneal wounds to heal.

### RGTA (MATRIX THERAPY)

- Replaces and mimicks heparan sulphate
- Preserves extracellular matrix (ECM) allowing normal tissue regeneration to take place.
- Animal discomfort reduced with improved quality ECM.

### INSTRUCTIONS: FOR OCULAR USE

- Use 1 drop every second day for small animals.
- Use 2 drops every second day for equine & bovine.
- After applying, take a 2-3 minute interval between the use of Clerapliq and other topical eye products.
- Store between 4-25 degrees Celsius.

For further information,  
video on mechanism  
of action and access  
to information for use,  
please scan QR code.





# A great new solution to an age-old problem...

## Rear end instability

### GIVE NEW LIFE TO OLD PAWS

Senior and disabled dogs lose traction and the ability to navigate smooth surfaces like they used to.

PawFriction pad coating doesn't require individual sizing and is easy to apply. In just a few minutes, you can give your dog a new lease on life and the ability to navigate smooth surfaces like they used to.

### PAWFRICITION IS MOST EFFECTIVE FOR DOGS WITH:


- Arthritis and old age.
- Neurological conditions: degenerative myelopathy, disc disease, lumbosacral instability, brain tumours, Wobbler's
- Orthopaedic issues: fractures, ACL injuries, patella luxation
- Post-operative recovery: back surgery, knee surgery (TPLO, TTA, lateral suture), amputations
- Confidence issues: reluctant to walk on smooth floors
- Rehabilitation or physical therapy needs
- Other conditions that can make walking on smooth floors difficult

### PAWFRICITION KIT CONTAINS

- 4x applicators of Methyl cyanoacrylate viscous glue
- 2x cannisters of fine rubber particles (non-toxic)

### DURABILITY

- One application lasts 1 week on a larger dog (4 weeks to a pack)
- For small dogs a pack will last around 15 weeks.




Giving new life to old paws Simple to apply, proven to work			
CODE	DESCRIPTION	QUANTITY	PRICE
PAW001	PawFriction Kit (sufficient for 16 paws)	1	\$65

# Vital Signs Monitor

## Vetcorder

Vital contact for every companion			
CODE	DESCRIPTION	QUANTITY	PRICE
VC-48525	Vital Signs Monitor	1	\$1750

### Reach New Levels of Care, Connection & Convenience




## FREE MOBILE/TABLET APP WHICH PROVIDES COMPLETE CONNECTION WITH VET CORDER, MIRRORING THE MEASUREMENT IN REAL TIME.

### MONITOR ANYWHERE, ANYTIME ON YOUR PHONE OR TABLET

Expand your practice beyond the four walls of your hospital. Use the Vetcorder in surgery, in your mobile practice or send it home with your client and receive regular vital signs updates.


### DO MORE TO MONITOR EVERY PET'S RECOVERY

The Vetcorder offers an alternative to outdated, high cost and bulky monitors. As a portable smart monitor, it stays with your patient from admission through discharge.



Battery Life
Over 12 hours of continuous monitoring
3 months standby

SpO2	
SPO2 RANGE	70% to 100%
SpO2 Accuracy (Arms)	80-100% 2%, 70-79%: 3%
PR Range	30 to 360 bpm
PR Accuracy	±2 bpm or ±2%, whichever is greater





# Patient Warming

## Hot Dog®

### PREVENT AND TREAT HYPOTHERMIA IN SURGICAL PATIENTS

CODE	DESCRIPTION	QUANTITY	PRICE
HDSP-17	HotDog® Single Port Controller	1	\$3,000
HDDP-17	HotDog® Dual Port Controller	1	\$4,300
HDBXS	HotDog® Blankets X-Small 10' x 12'	1	\$549
HDBS	HotDog® Blankets Small 10' x 23'	1	\$585
HDBM	HotDog® Blankets Medium 16' x 22'	1	\$625
HDBL	HotDog® Blankets Large 22' x 31'	1	\$697
HDBXL	HotDog® Blankets Extra Large 27' x 47'	1	\$742
HDBC	Blanket Cable	1	\$160

## Hot Dog® – No hot air... No water... Just safe, even warmth!

**HOTDOG® IS A THERMOSTATICALLY CONTROLLED, HIGHLY EFFECTIVE, INTRA-OPERATIVE WARMING SYSTEM.**

**BY NOT BLOWING FORCED, WARMED AIR INTO THE O.R., ANY RISKS ASSOCIATED WITH THIS PRACTICE ARE AVOIDED.**

- No forced warm air in the O.R.
- No ongoing consumables
- Durable, just wipe down the blankets after use
- Thermostatically controlled
- Highly effective warming
- Safe



### Comparison with other elements

FEATURES	HOT DOG®	HOT AIR (BAIR HUGGER)	HOT WATER
Effective warming	✓	✓	
Even heat distribution	✓		
Low ongoing costs	✓		
Durable waterproof blankets	✓		
Easy to clean	✓		
Won't dry EKG leads or mucous membranes	✓		✓
Maintenance-free controller	✓		
Quiet operation	✓		✓
Fast warm-up time	✓	✓	
Small and light for ease of use	✓		
No blowing air on surgical site	✓		✓
Easily fits through cage slats	✓		✓
Won't incubate bacteria	✓		
Environmentally friendly	✓		✓
Energy efficient	100 watts	800 watts	Up to 800 watts




SCAN ME




# Needle Chart

## Needle Shapes

1/2 Circle




3/8 Circle




Straight Cutting




1/2 Circle Conventional Cutting /  
Straight Shank (Ski)




## Needle Types



Round bodied taper point needle ●













Reverse cutting needle ▽




Conventional cutting needle △

# Circle Taper Point

## 1/2 Circle Taper Point

Needle	Size	SilverGlide™	Ethicon™ (similar)	Syneture™	B. Braun™ (similar)
	• 11mm	T11	-		
	• 13mm	T22	TF	CV-22	HR13
	• 16mm	T16	-		
	• 18mm	T28	RB-1	CV-23	HR18
	• 22mm	T4	SH-1	CV-25	HR22
	• 26.5mm	T5	SH / CT-1	V20 / GS22	HR26
	• 36mm	T8	MH / V-26	V-26	HR36
	• 37mm	T37	CT-1	-	-
	• 40mm	T40	-	-	-
	• 65mm	T39	TP-1	-	-

## 3/8 Circle Taper Point


Needle	Size	SilverGlide™	Ethicon™ (similar)	Syneture™	B. Braun™ (similar)
	• 90mm	T90	-	-	-

# Circle Reverse Cutting

## 1/2 Circle Reverse Cutting ▽

Needle	Size	SilverGlide™	Ethicon™ (similar)	Syneture™	B. Braun™ (similar)
	▽ 13mm	C13	-	-	-
	▽ 26mm	C21	CP-2	GS-10	HS26
	▽ 37mm	C37	CP-1	GS-11	HS37
	▽ 40mm	C14	CP	GS-12	HS40

## 3/8 Circle Reverse Cutting ▽

Needle	Size	SilverGlide™	Ethicon™ (similar)	Syneture™	B. Braun™ (similar)
	▽ 11mm	C16	C-2	C-1	DS12
	▽ 12mm	C17	C-2	C-1	DS12
	▽ 16mm	C22	FS-3	C-12	DS16
	▽ 18mm	C6	FS-2	C-13	DS18
	▽ 24mm	C7	FS-1	C-14	DS24
	▽ 25mm	C25	-	-	-
	▽ 26mm	C26	FS	-	-
	▽ 30mm	C10	FSL	C-16	DS30
	▽ 40mm	C14	FSLX	C-17	
	▽ 50mm	C50	-	-	-
	▽ 77mm	C77	LR	-	-













# Straight Cutting

## Straight Cutting ▽

Needle	Size	SilverGlide™	Ethicon™ (similar)	Syneture™	B. Braun™ (similar)
 Straight Cutting	▽ 60mm	CE60	KS	SC-2	GS60
 Straight Cutting	▽ 65mm	CE65	KS	SC-2	GS60
 1/2 Circle Conventional Cutting / Straight Shank (Ski) △	△ 75mm	SKI75	-	-	-



# Suture Material Absorption Profile

NON-ABSORBABLE		MATERIAL DESCRIPTION	QUALITIES		SIMILAR TO
NYLON		Black monofilament	High strength retention		Ethilon™ * Monosor™ §/ Dermalon™ § Dafilon™ †
FLUOROMID™		Fluorescent Yellow polyamide multifilament	Easy handling, ultra high visibility Fluorescent yellow		
WHITE FLUOROMID™		White polyamide multifilament	Easy handling		
SHORTY™ POLYVINYLIDENE FLUORIDE YELLOW		Bright Yellow Monofilament	High visibility Exceptional Knot Security, knots flatten and lock against themselves, Unwetted by blood, inert		
BLUE POLYPROPYLENE		Dark blue monofilament	Unwetted by blood, inert, knots flatten and lock against themselves for improved knot security		Prolene™ * Surgipro™ §/ Surgilene™ § Premilene™ †
FLUORESCENT PINK POLYPROPYLENE		Fluorescent Pink monofilament	Ultra high visibility Fluorescent pink Unwetted by blood, inert, knots flatten and lock against themselves for improved knot security		
POLYESTER		Braided polyester green silicone coated	High strength braid		Ticron § Ethibond *
ABSORBABLE		MATERIAL DESCRIPTION	BSR PROFILE	MASS ABSORPTION	SIMILAR TO
GLYCOLIDE & ECAP		Monofilament purple	In vivo strength retention 60% at 7 days 25% at 14 days	Essentially complete between 2-4 months	Monocryl™ * Caprosyn™ §/ Biosyn™ § Monosyn™ †
POLYGLYCOLIC ACID		Braided violet	In vivo strength retention 65% at 14 days 20% at 21 days	Essentially complete between 3-6 months	Vicryl™ * Safil™ † Polsorb §
POLYGLYCOLIC ACID FAST ABSORBING		Braided undyed	In vivo strength retention 48% at 7 days 0% at 14 days	Essentially complete at 6 weeks	Vicryl Rapide™ * Safil Quick™ †
POLYGLACTIN		Braided violet	In vivo strength retention 74% at 14 days 49% at 21 days 24% at 28 days	Essentially complete at 9 weeks	Vicryl™
POLYDIOXANONE		Monofilament purple	In vivo strength retention 80-90% at 14 days 60-80% at 28 days 50-80% at 42 days	Essentially complete between 6-9 months	PDS11™ * MonoPlus™ † Maxon §

\* Trademark of Ethicon   § Trademark of Covidien   † Trademark of B.Braun





For sample requests or to place an order call  
Carolyn on **1300 365 241** or visit our website

***[silverglide.com.au](http://silverglide.com.au)***

**Distributed by:**